

**Category:** Biochemistry

**Title:** Estimation of Glucose on Mini CHEM 100 using serum samples

**SOP No.:** 10/03

**Date first effective:** 1<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

**Category:** Biochemistry

**Title:** Estimation of Glucose in serum samples of patients and normal, healthy volunteers on Mini CHEM 100

**SOP No:** 10/03

**Total pages:** 4

**Date first effective:** 1<sup>st</sup> January 2025

**Next Review date:** 31<sup>st</sup> December 2025

**Author:** Vaishali Hadikar  
Laboratory Technician

Signature with date:

*Hadikar*  
31/Dec/24

Department of Clinical Pharmacology  
Seth GS Medical College & KEM Hospital,  
Parel, Mumbai - 400 012.

**Lab In charge:** Dr. Sheetal Kudtarkar  
Project Scientist

Signature with date:

*Sheetalkar*  
31/Dec/24

**Reviewer:** Dr. Mahesh Belhekar  
Associate Professor

Signature with date:

*Belhekar*  
31/Dec/24

**Mahesh N. Belhekar**  
Associate Professor  
Department of Clinical Pharmacology  
New MS Building, First Floor,  
Seth GS Medical College and KEM Hospital  
Acharya Donde Marg, Parel,  
Mumbai - 400 012, India

**Approved by:** Dr. Nithya Gogtay  
Professor and Head

Signature

Confidential

*Nithya*  
31/12/24  
**Dr. Nithya Gogtay**  
Professor & Head  
Department of Clinical Pharmacology  
1<sup>st</sup> Floor, MS Building,  
Seth GS Medical College & KEM Hospital  
Parel, Mumbai - 400 012.

# CURRENT

Category: Biochemistry

Title: Estimation of Glucose on Mini CHEM 100 using serum samples

SOP No.: 10/03

Date first effective: 1<sup>st</sup> January 2025

Review date: 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

## Table of

## contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	References to other applicable SOPs	3
5	Detailed Instructions	3

### 1. Purpose

**Category:** Biochemistry

**Title:** Estimation of Glucose on Mini CHEM 100 using serum samples

**SOP No.:** 10/03

**Date first effective:** 1<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

The purpose of this Standard Operating Procedure (SOP) is to outline the procedure for estimation of glucose on Mini CHEM 100 in the Biochemistry lab of Department of Clinical Pharmacology.

## 2. Scope

This SOP covers the procedure of estimation of glucose on Mini CHEM 100 in serum.

## 3. Responsibility

Lab technician, Lab attendant, or any other appropriately qualified staff in the team, designated by the Head of Department, will be responsible for analysis.

## 4. Reference

- Departmental SOP no 10/05: Blood collection.
- Departmental SOP no 24/04: Waste management.
- Biochemistry labs SOP no 10/03: Mini CHEM 100 instrument.
- Kit insert

## 5. Detailed Instructions

1. The whole blood sample is collected in Sodium fluoride bulb as per departmental SOP no. 10/04 and is processed as per biochemistry lab SOP no 01/04.
2. Remove the kit kept in the **Refrigerator** located in the main Biochemistry laboratory of Department of Clinical Pharmacology, M.S. Building, 1<sup>st</sup> Floor.

**Use kit:-** SPIN REACT

**KIT Expiry:** 2 years after the manufacturing date mentioned on kit



# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Glucose on Mini CHEM 100 using serum samples

**SOP No.:** 10/03

**Date first effective:** 1<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

Manufactured by: - SPIN REACT

**Reagent 1: Glucose (Reagent is ready to use)**

**Standard Solution: 100mg/dl**

- The unopened reagent remains stable till expiry date at 2-8°C, when protected from contamination.
- 3. As per the instruction displayed on the machine, aspirate distilled water for washing purpose taken from Milli Q water purifier .
- 4. Press **GLU** displayed on the machine.
- 5. When Glucose parameters are displayed, water is aspirated, followed by reagent blank and sample in glass tubes. (the tubes uses for Mini CHEM 100 are kept in a labeled box on left side of the machine).
- 6. Dilution is prepared by adding **10µl** of distilled water in **1000µl** of Glucose reagent in new glass tube.
- 7. Add **10µl** of standard in **1000µl** of Glucose reagent in new glass tube.
- 8. **Add 10µl** of sample in **1000µl** of Glucose reagent in new glass tube.
- 9. For more than one sample, **add 10µl** of sample in **1000µl** of Glucose reagent in new glass tube.
- 10. Mix well and incubate for 15 minute at 37°C in an incubator in biochemistry lab.
- 11. Aspirate the standard followed by the samples( Glucose standard: 100 mg/dl)
- 12. Calculation is done automatically by the machine. Glucose activity [**mg/dl**]
- 13. **Linearity:** Upto 500mg/dl. For higher values it is recommended to dilute the samples with normal saline and repeat the assay. Multiply the results with the dilution factor.