

Category: Biochemistry

Title: Estimation of Calcium in blood serum samples on Mini CHEM 100

SOP No.: 11/03

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

Category: Biochemistry

Title: Estimation of calcium in blood serum samples on Mini CHEM 100

SOP No: 11/03

Total pages: 5

Date first effective: 01st January 2025

Next Review date: 31st December 2025

Author: Vaishali Hadikar
Laboratory Technician

Signature with date:

Hadikar
31 Dec 2024

Department of Clinical Pharmacology
Seth GS Medical College & KEM Hospital,
Parel, Mumbai - 400 012.

Lab In charge: Dr. Sheetal Kudtarkar
Project Scientist

Signature with date:

Sheetal Kudtarkar
31 Dec 24

Reviewer: Dr. Mahesh Belhekar
Associate Professor

Signature with date

Belhekar
31 Dec 2024

Dr. Mahesh N. Belhekar
Associate Professor
Department of Clinical Pharmacology
New MS Building, First Floor,
Seth GS Medical College and KEM Hospital
Acharya Dadas Marg, Parel,
Mumbai - 400 012, India

Approved by: Dr. Nithya Gogtay
Professor and Head

Signature with date:

Nithya Gogtay
31.12.24

Dr. Nithya Gogtay
Professor & Head
Department of Clinical Pharmacology
1st Floor, MS Building,
Seth GS Medical College & KEM Hospital,
Parel, Mumbai - 400 012.

CURRENT

Category: Biochemistry

Title: Estimation of Calcium in blood serum samples on Mini CHEM 100

SOP No.: 11/03

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

Table of contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	References to other applicable SOPs	3
5	Detailed Instructions	3

Category: Biochemistry

Title: Estimation of Calcium in blood serum samples on Mini CHEM 100

SOP No.: 11/03

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

1. Purpose

The purpose of this Standard Operating Procedure (SOP) is to outline the procedure for estimation of calcium on Mini CHEM 100 in the Biochemistry lab of Department of Clinical Pharmacology

2. Scope

This SOP covers the procedure of estimation of calcium on Mini CHEM 100 in blood serum.

3. Responsibility

Lab technician, Lab attendant, or any other appropriately qualified staff in the team, designated by the Head of Department, will be responsible for analysis.

4. Reference

- Departmental SOP no 10/05: Blood collection.
- Departmental SOP no 24/04: Waste management.
- Biochemistry labs SOP no 10/03: Mini CHEM 100 instrument.
- Kit insert

5. Detailed Instructions

1. The whole blood sample is collected in plain bulb as per departmental SOP no. 10/04 and is processed as per biochemistry lab SOP no 01/04.
2. Remove the kit kept in the **Refrigerator** located in the main Biochemistry laboratory of Department of Clinical Pharmacology, M.S. Building, 1st Floor.

Use kit: SPINREACT Calcium KIT,

CURRENT

Category: Biochemistry

Title: Estimation of Calcium in blood serum samples on Mini CHEM 100

SOP No.: 11/03

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

KIT Expiry: 2 years after the manufacturing date mentioned on kit

Manufactured by: - Spinreact

3. **Reagent 1: Calcio Arsenzo (Reagent is ready to use)**

Working Reagent Stability:-

- The unopened reagent is stable till expiry date when stored at (2-8°C), when working solution become turbid discard it.
4. As per the instruction displayed on the machine, aspirate distilled water for washing purpose taken from Milli Q water purifier of Post PCR Lab/.
 5. Press **CAL** displayed on the machine.
 6. When calcium parameters are displayed, water is aspirated, followed by reagent blank, standard and sample in glass tubes (the tubes uses for Mini CHEM 100 are kept in a labeled box on left side of the machine).
 7. Dilution is prepared for reagent blank by adding **10µl** of distilled water in **500µl** of calcium arsenazo reagent in new glass tube.
 8. Add **10µl** of standard in **500µl** of calcium arsenazo reagent in new glass tube.
 9. Add **10µl** of sample in **500µl** of calcium arsenazo reagent in new glass tube.
 10. For more than 1 sample, add **10µl** of sample in **500µl** of calcium arsenazo reagent in new glass tube.
 11. Mix well and read the absorbance against blank.
 12. Aspirate the standard followed by the samples (Calcium standard solution: 10.0 mg/dl)
 13. Calculation is done automatically by the machine. Calcium acid activity [**mg/dl**].
 14. **Linearity:** Upto 16 mg/dl. For higher values it is recommended to dilute the samples with normal saline and repeat the assay. Multiply the results with the dilution factor.
 15. The results are entered in form sent with the sample by the concerned doctor. The form is attached in Appendix 1 as in biochemistry lab SOP no 10/03. The results are also entered in biochemistry report book.
 16. Values are evaluated and signed by biochemist, /laboratory in charge and Study-coordinator and PI.

Category: Biochemistry

Title: Estimation of Calcium in blood serum samples on Mini CHEM 100

SOP No.: 11/03

Review date: 31st December 2025

Date first effective: 1st January 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

17. Internal Quality control is done once in a week as per biochemistry labs SOP no 05/01
(SOP for external and internal QCs).
18. Once the calcium parameter is completed for all samples, aspirate distilled water for washing purpose
19. After washing is completed, switch off the adaptor and then switch off the main switch.
20. Samples are discarded as per departmental SOP no. 24/04 (SOP of waste management).