

# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Alkaline Phosphatase (ALP) in blood serum samples on minichem ARK Diagnostics

**SOP No.:** Biochemistry 4/02

**Date first effective:** 01<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

**Category:** Biochemistry

**Title:** Estimation of Alkaline phosphatase (ALP) in blood serum samples on Minichem ARK Diagnostics

**SOP No:** Biochemistry04/02

**Total pages:** 4

**Date first effective:** 02<sup>nd</sup> January 2025

**Next Review date:** 31<sup>st</sup> December 2025

**Author:** Vaishali Hadikar  
Laboratory Technician

Signature with date:

*V Hadikar*  
31/Dec/2024

**Lab In charge:** Dr. Sheetal Kudtarkar  
Project Scientist

Signature with date:

*Sheetal Kudtarkar*  
31/Dec/24

**Reviewer:** Dr. Mahesh Belhekar  
Associate Professor  
Department of Clinical Pharmacology  
New MS Building, First Floor,  
Seth GS Medical College and KEM Hospital  
Acharya Donda Marg, Parel,  
Mumbai - 400 012. India

Signature with date

*M. Belhekar*  
31/Dec/2024

**Approved by:** Dr. Nithya Gogtay  
Professor and Head

Signature with date:

*Nithya Gogtay*  
31/12/24  
Dr. Nithya Gogtay  
Professor & Head  
Department of Clinical Pharmacology  
1<sup>st</sup> Floor, MS Building,  
Seth GS Medical College & KEM Hospital  
Parel, Mumbai - 400 012.

# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Alkaline Phosphatase (ALP) in blood serum samples on minichem ARK Diagnostics

**SOP No.:** Biochemistry 4/02

**Date first effective:** 01<sup>st</sup>January2025

**Review date:** 31stDecember 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

## Table of contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	References to other applicable SOPs	3
5	Detailed Instructions	3

**Category:** Biochemistry

**Title:** Estimation of Alkaline Phosphatase (ALP) in blood serum samples on minichem ARK Diagnostics

**SOP No.:** Biochemistry 4/02

**Date first effective:** 01<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

## 1. Purpose

The purpose of this Standard Operating Procedure (SOP) is to outline the procedure for estimation of Alkaline Phosphatase (ALP) on minichem ARK Diagnostics in the Biochemistry lab of Department of Clinical Pharmacology

## 2. Scope

This SOP covers the procedure of estimation of Alkaline Phosphatase (ALP) on minichem ARK Diagnostics.

## 3. Responsibility

Lab technician, Lab attendant, or any other appropriately qualified staff in the team, designated by the Head of Department, will be responsible for analysis.

## 4. Reference

- Departmental SOP no 10/05: Blood collection.
- Departmental SOP no 24/04: Waste management.
- Biochemistry labs SOP no 01/01: Analysis on minichem ARK diagnostics.
- Kit inserts

## 5. Detailed instructions

1. The whole blood sample is collected and processed as per SOP no. 01/03.
2. Remove the Spin react reagent kit kept in the **Refrigerator** [between 2-8 degree Celsius] located in the main Biochemistry laboratory of Department of Clinical Pharmacology, M.S. Building, 1<sup>st</sup> Floor. Verify the expiry date before using the kit.  
**Use kit:** Spin react ALP Kit,
3. As per SOP no. 10/03 for Biochemistry analysis on minichem ARK Diagnostics, select the parameter 'ALP' and set the machine for measurement of the sample.
4. As per the reagent kit insert, to 10 µl of the sample add 480 µl of reagent 1 and 120 µl of reagent 2.

# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Alkaline Phosphatase (ALP) in blood serum samples on minichem ARK Diagnostics

**SOP No.:** Biochemistry 4/02

**Date first effective:** 01<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

5. After addition of the reagents, immediately aspirate the test solution through the sipping tube by pressing the button behind the tube.
6. The machine will give a beep sound and will continue with the measurement of the concentration for a period of 120 seconds.
7. After the delay time of 120 seconds, the machine will beep and will display the absorbance followed by the reading.
8. Values are entered in the report, evaluated and signed by the biochemist /laboratory in-charge, Study-coordinator and PI.
9. Once the ALP parameter is completed for all samples, aspirate distilled water for washing purpose.
10. After washing is completed, switch off the adaptor and then switch off the main switch.
11. Samples are discarded as per departmental SOP no. 24/02 (SOP of waste management).