

# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Creatinine on minichem ARK diagnostic using serum samples

**SOP No.:** Biochemistry05/01

**Date first effective:** 01<sup>st</sup>January2025

**Review date:** 31<sup>st</sup>December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

**Category:** Biochemistry

**Title:** Estimation of creatinine in serum samples of patients and normal, healthy  
Volunteers on minichem ARK diagnostic

**SOP No:** Biochemistry 05/01

**Total pages:** 4

**Date first effective:** 01<sup>st</sup>January2025

**Next Review date:** 31<sup>st</sup>December 2025

**Author:** Vaishali Hadikar  
Laboratory Technician

Signature with date:

*Hadikar*  
31 Dec 2024

Department of Clinical Pharmacology  
Seth GS Medical College & KEM Hospital,  
Parel, Mumbai - 400 012.

**Lab In charge:** Dr. Sheetal Kudtarkar  
Project Scientist

Signature with date:

*Sheetal Kudtarkar*  
31/Dec/24

**Reviewer:** Dr. Mahesh Belhekar  
Associate Professor

Signature with date

*Belhekar*  
31/Dec/24

**Mahesh N. Belhekar**  
Associate Professor  
Department of Clinical Pharmacology  
New MS Building, First Floor,  
Seth GS Medical College and KEM Hospital  
Acharya Dande Marg, Parel,  
Mumbai - 400 012, India

**Approved by:** Dr. Nithya Gogtay  
Professor and Head

Signature with date:

*Nithya Gogtay*  
31/12/24

**Dr. Nithya Gogtay**  
Professor & Head  
Department of Clinical Pharmacology  
1<sup>st</sup> Floor, MS Building,  
Seth GS Medical College & KEM Hospital  
Parel, Mumbai - 400 012.

# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Creatinine on minichem ARK diagnostic using serum samples

**SOP No.:** Biochemistry05/01

**Date first effective:** 01<sup>st</sup>January2025

**Review date:** 31<sup>st</sup>December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

## Table of contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	References to other applicable SOPs	3
5	Detailed Instructions	3

**Category:** Biochemistry

**Title:** Estimation of Creatinine on minichem ARK diagnostic using serum samples

**SOP No.:** Biochemistry05/01

**Review date:** 31<sup>st</sup>December 2025

**Date first effective:** 01<sup>st</sup>January2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

## 1. Purpose

The purpose of this Standard Operating Procedure (SOP) is to outline the procedure for estimation of creatinine on minichem ARK Diagnostics semi automated machine in the Biochemistry lab of Department of Clinical Pharmacology.

## 2. Scope

This SOP covers the procedure of estimation of creatinine on minichem ARK Diagnostics semi automated machine using serum.

## 3. Responsibility

Lab technician, Lab attendant, or any other appropriately qualified staff in the team, designated by the Head of Department, will be responsible for analysis.

## 4. Reference

- Departmental SOP no 10/05: Blood collection.
- Departmental SOP no 24/04: Waste management.
- Biochemistry labs SOP no 1/01: Analysis using minichem ARK diagnostic instrument.
- Kit inserts

## 5. Detailed Instructions

1. The whole blood sample is collected and processed as per SOP no. 10/05.
2. Remove the QCA reagent kit kept in the **Refrigerator**[between 2-8 degree Celsius] located in the main Biochemistry laboratory of Department of Clinical Pharmacology, M.S. Building, 1<sup>st</sup> Floor. Verify the expiry date before using the kit.
3. **Use kit:** Spin react Jaffe's kinetic method Creatinine Kit,

# CURRENT

**Category:** Biochemistry

**Title:** Estimation of Creatinine on minichem ARK diagnostic using serum samples

**SOP No.:** Biochemistry05/01

**Date first effective:** 01<sup>st</sup>January2025

**Review date:** 31<sup>st</sup>December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012

4. As per SOP no. 10/03 for Biochemistry analysis on minichem ARK Diagnostics, select the parameter 'CREATININE' and set the machine for measurement of the sample.
5. As per the reagent kit insert, to 50 µl of the sample add 250µl of reagent 1 and 250µl of reagent 2.
6. After addition of the reagents, immediately aspirate the test solution through the sipping tube by pressing the button behind the tube.
7. The machine will give a beep sound and will continue with the measurement of the concentration for a period of 120 seconds.
8. After the delay time of 120 seconds, the machine will beep and will display the absorbance followed by the reading.
9. Values are entered in the report, evaluated and signed by the biochemist/laboratory in charge and Study-coordinator and PI.
10. Once the Creatinine parameter is completed for all samples, aspirate distilled water for washing purpose.
11. After washing is completed, switch off the adaptor and then switch off the main switch.
12. Samples are discarded as per departmental SOP no. 24/02 (SOP of waste management).