

CURRENT

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.

Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Category: General SOP for Extraction of genomic DNA from whole blood.

Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Parel, Mumbai.

SOP No.: Pre PCR-Method- 06

Total pages: 07

Date first effective: 1st January 2025

Review date: 31st December 2025

Version: 06

Author: Dr. Sheetal Kudtarkar
Project Associate

Signature with date:

Sheetal
30/12/24

Reviewer 1: Dr. Bhaskar Krishnamurthy
Assistant Professor

Signature with date:

Bhaskar JK
30/DEC/2024

Dr. Bhaskar Krishnamurthy
Assistant Professor,

Department of Clinical Pharmacology,
Seth GSMC and KEMH, Mumbai -400 012.

Reviewer 2: Dr. Mahesh Belhekar
Associate Professor

Signature with date:

Belk
30/Dec/2024
Dr. Mahesh N. Belhekar
Associate Professor
Department of Clinical Pharmacology
New MS Building, First Floor,
Seth GS Medical College and KEM Hospital
Acharya Donde Marg, Parel,
Mumbai - 400 012, India

Approved by: Dr. Nithya Gogtay
Professor and Head

Signature with date:

Dr. Nithya Gogtay

Professor & Head
Department of Clinical Pharmacology
1st Floor, MS Building,
Seth GS Medical College & KEM Hospital
Parel, Mumbai - 400 012.

Nithya
30/12/24

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.

Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Table of contents

No.	Contents	Page No.
1.	Purpose	3
2.	Scope	3
3.	Responsibilities	3
4.	References to other applicable SOPs	3
5.	Detailed Instructions	4
6.	Picture	6 -7

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.

Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Purpose: This standard operating procedure (SOP) describes the method for extraction of genomic DNA from whole blood by kit method

Scope: This SOP is limited for extraction of genomic DNA from whole blood by kit method

Responsibilities:

Primary Responsibility: Divya Bhare
Lab Technician

Divya Bhare
30/Dec/2024

Secondary Responsibility: Dr. Sheetal Kudtarkar
Project Scientist

Sheetal Kudtarkar
30/12/24

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.

Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Detailed Instructions:

Important: This procedure should be carried out wearing sterile gloves

Reagents preparation:

- Proteinase K: Add 1.34 mL double distilled water to lyophilized Proteinase K powder.
- Wash buffer 1: Add 8mL absolute ethanol to concentrated wash buffer1 to make working wash buffer.
- Wash buffer: Add 40mL absolute ethanol to concentrated wash buffer to make working wash buffer.

Procedure:

1. Pipette 20 μ L proteinase K into the bottom of a 1.5 mL micro -centrifuge tube, to it add 200 μ L blood/ body fluid sample or buffy coat from 1ml blood (equilibrated to room temperature).
2. Add 200 μ L buffer FABG to the samples and vortex the mixture vigorously (30 sec). Incubate samples at 70^oC for 10-15 mins.
3. Add 200 μ L ethanol (96-100%) to each sample and vortex again.
4. Pipette the mixture from step 3 into the DNA Sure Blood Mini Column placed in a collection tube. Centrifuge 1 min at 11000 \times g (see below). Discard collection tube with flow through.
5. Place the FABG Mini Column into fresh collection tube (2mL) and add 400 μ L W1 Buffer. Centrifuge 1 min at 11000 \times g. Discard the collection tube with flow through.
6. Place the FABG Mini Column into fresh collection tube (2mL) and add 700 μ L WASH buffer. Centrifuge 1 min at 11000 \times g. Discard flow through and reuse collection tube.
7. Place the FABG Mini Column back into the collection tube and centrifuge for 1 min at 11000 \times g
8. Place the FABG Mini Column in a 1.5mL micro-centrifuge tube and add 100 μ L preheated elution buffer (70^oC). Dispense buffer directly onto silica membrane. Incubate at room temperature for 1min. Centrifuge for 1min at 11000 \times g.
9. Flow through contains the DNA. Store at 4^oC.

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.

Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Note: The Relative Centrifugal Force (RCF) is expressed in units of gravity (time's gravity or $\times g$). Many microcentrifuges only have settings for speed (revolutions per minute, RPM), not relative centrifugal force. Consequently, a formula for conversion is required to ensure that the appropriate setting is used in an experiment. The relationship between RPM and RCF is as follows:

$$g = (1.118 \times 10^{-5}) r S^2$$

Where 'g' is the relative centrifugal force, 'r' is the radius of the rotor in centimeters, and 'S' is the speed of the centrifuge in revolutions per minute.

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.


Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.



DNA extraction room



Reagents used

Category: General SOP for Extraction of genomic DNA from whole blood in Genetics laboratory.


Title: Extraction of genomic DNA from whole blood by kit method in the Department of Clinical Pharmacology, K.E.M Hospital, Mumbai.

SOP No.: Pre PCR-Method- 06

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st floor, New MS building Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.



Cooling centrifuge



Column & Collection tube used for DNA extraction