

CURRENT

Category : Patient Care
Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM010/01

Date first effective: 1st January 2025 Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Category : Patient Care

Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No.: TDM010

Total pages: 6

Date first effective: 1st January 2025

Next Review date: 31st December 2025

Version: 01

Author:

Swati More
30/Dec/24
Swati More
Lab. Analyst

Dr. Bhaskar Krishnamurthy
Lab. In-charge

Signature with date:

BK
30/DEC/2024

Reviewer:

Belk
30/Dec/2024
Dr. Mahesh Belhekar
Associate Professor

Dr. Bhaskar Krishnamurthy
Assistant Professor,
Department of Clinical Pharmacology,
Seth GSMC and KEMH, Mumbai -400 012.

Signature with date:

Dr. Mahesh M. Belhekar
Associate Professor
Department of Clinical Pharmacology
New MS Building, First Floor,
Seth GS Medical College and KEM Hospital
Acharya Donde Marg, Parel,
Mumbai - 400 012. India

Approved by:

Dr. Nithya Gogtay
Professor and Head

Signature with date:

N
30/12/24
Dr. Nithya Gogtay
Professor & Head
Department of Clinical Pharmacology
1st Floor, MS Building,
Seth GS Medical College & KEM Hospital,
Parel, Mumbai - 400 012.

Confidential

Page 1 of 6



Category : Patient Care
Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM010/01

Date first effective: 1st January 2025 **Review date:** 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
 Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

1. Table of Contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	Applicable rules/guidelines	3
5	Reference to other applicable SOPs	3
6	Detailed instructions	4
7	Abbreviations	6

Category : Patient Care
Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM010/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

1. Purpose:

This SOP describes the technique for qualitative and quantitative estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

2. Scope:

This SOP is limited to estimation of Vancomycin in $\mu\text{g/ml}$ in human serum by High Performance Liquid Chromatography (HPLC).

3. Responsibilities:

The Head of the department is responsible for the medical care and welfare of all patients pertaining to TDM of Vancomycin. The task of performing estimation of Vancomycin will be delegated to trained personnel who will perform this function.

4. Applicable rules, regulations and guidelines

- ICMR Good Clinical Laboratory Practices Guidelines 2021 (<http://icmr.nic.in/guidelines/GCLP.pdf>)

5. Reference to other applicable SOPs

- SOP No.24/02: Biomedical waste management.
- SOP No. TDM01/02: Collection and separation blood plasma for TDM
- SOP No. TDM05/02: Operation of High-Performance Liquid Chromatography

Category : Patient Care
Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM010/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

6. Detailed instructions

1. Preparation of standard and calibrator

a. Preparation of Stock Standard (Vancomycin 1mg/mL) : (10mg)

10 mg of pure powder of Vancomycin + 10 mL of distilled water

b. Preparation of stock Internal standard (IS): [Caffeine (CAF)]

10mg of pure powder (Caffeine) IS + 10mL of Distilled water-stock IS

c. Preparation of Working standard (IS): [Caffeine (CAF)]

Working standard I	1mL of stock IS	9mL of Distilled water
Working standard II	5 ml of working standard (I)	5 mL of Distilled water
Working standard III	5 ml of working standard (II)	5 mL of Distilled water.

d. Preparation of Mobile Phase buffer (50mM):

Weigh 5.7515g of Ammonium dihydrogen phosphate ($\text{NH}_4\text{H}_2\text{PO}_4$) and dissolve in 1000 mL of distilled water

2. Preparation of serum standards:

Concentration	Plasma	Vancomycin	Final volume
20 $\mu\text{g/mL}$	980 μL	20 μL of 1mg/ml Stock	1000 μL
10 $\mu\text{g/mL}$	500 μL	500 μL of 10 $\mu\text{g/ml}$ Stock	1000 μL
5 $\mu\text{g/mL}$	500 μL	500 μL of 20 $\mu\text{g/ml}$ Stock	500 μL

Category : Patient Care
Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM000/01

Date first effective: 1st January 2025 **Review date:** 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

3. Extraction procedure

500 μ L of serum (blank/standards/quality control/patients' sample) in
2mL Eppendorf tube.



Add 30 μ L of internal standard (Working standard III) to it, Vortex mix
thoroughly for 30 seconds.



Add 500 μ L of Isopropanol: Acetonitrile (1:1) diluent and Vortex mix
thoroughly for 2 min.



Keep the Eppendorf tubes in ice tub for 10 minutes to facilitate protein
precipitation.



Remove from ice tub after 10 minutes and shake 100 times manually and again
vortex for 2 minutes.



Centrifuge the suspension at 15000 rpm for 15mins at 4°C.



Transfer the Supernatant with Pasteur pipette to another 7mL glass test tubes.



Add 3.5mL methylene chloride to the supernatant and vortex for 2 minutes.

Again, centrifuge at 2500 rpm for 20 minutes



Remove the 100 μ L of supernatant and add transfer into HPLC vials.

Inject into HPLC for analysis.

Category : Patient Care
Title : Procedure for estimation of Vancomycin in human serum by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM010/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

4. Preparation of Mobile phase

- a. Take Ammonium dihydrogen phosphate ($\text{NH}_4\text{H}_2\text{PO}_4$) buffer (50mM): Acetonitrile; in the ratio of 90:10 and adjust the PH to 2.2 with orthophosphoric acid and 100 μL formic acid for 1000mL.
- b. Filter the mobile phase through 0.22-micron filter and sonicate for 15 minutes.

5. HPLC Conditions

- a. Injecting volume: 25 μL
- b. Flow rate: 0.90 mL/min.
- c. Wavelength: 205 nm (UV detector)
- d. Run Time: 20.00 min (approximately)
- e. Retention times for VAN- 10.0-11.0 min, I.S – 11.0– 12.0 min approximately.
- f. HPLC Column: Supelco C18 column (250 \times 4.6mm, 5 μL).

6. Abbreviations:

- a. HPLC = High Performance Liquid Chromatography
- b. I.S. = Internal Standard
- c. VAN = Vancomycin
- d. Std = Standard
- e. CAF = Caffeine

Confidential

Page 6 of 6

