

# CURRENT

Category : Patient Care  
Title : Procedure for estimation of Isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).  
SOP No. : TDM12/01  
Date first effective: 1<sup>st</sup> January 2025 Review date: 31<sup>st</sup> December 2025  
Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Category : Patient Care

Title : Procedure for estimation of Isoniazid in human plasma by  
Liquid Chromatography Mass Spectrometry (LCMS).

SOP No.: 12/01

Total pages: 10

Date first effective: 1<sup>st</sup> January 2025

Next Review date: 31<sup>st</sup> December 2025

Version: 01

Author:

*Swati More*  
30/Dec/24  
Swati More  
Lab. Analyst

Dr. Bhaskar Krishnamurthy  
Lab In-charge

Signature with date:

*BK*  
30/Dec/2024

Dr. Bhaskar Krishnamurthy  
Assistant Professor,  
Department of Clinical Pharmacology,  
Seth GSMC and KEMH, Mumbai -400 012.

Reviewer:

*Belh*  
30/Dec/2024  
Dr. Mahesh Belhekar  
Associate Professor

Signature with date:

Dr. Mahesh N. Belhekar  
Associate Professor  
Department of Clinical Pharmacology  
New MS Building, First Floor,  
Seth GS Medical College and KEM Hospital  
Acharya Donde Marg, Parel,  
Mumbai - 400 012. India

Approved by:

Dr. Nithya Gogtay  
Professor and Head

Signature with date

*N*  
30/12/24  
Dr. Nithya Gogtay  
Professor & Head  
Department of Clinical Pharmacology  
1<sup>st</sup> Floor, MS Building,  
Seth GS Medical College & KEM Hospital,  
Parel, Mumbai - 400 012.

Category : Patient Care  
 Title : Procedure for estimation of isoniazid in Human Plasma  
 by Liquid Chromatography Mass Spectrometry (LCMS).  
 SOP No. : TDM12/01  
 Date first effective: 1<sup>st</sup> January 2025      Review date: 31<sup>st</sup> December 2025  
 Department of Clinical Pharmacology, 1st Floor, New MS Building,  
 Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

### 1. Table of Contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	Applicable rules/guidelines	3
5	Reference to other applicable SOPs	3
6	Detailed instructions	4
7	Abbreviations	10

Category : Patient Care  
Title : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).  
SOP No. : TDM12/01  
Date first effective: 1<sup>st</sup> January 2025 Review date: 31<sup>st</sup> December 2025  
Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

1. Purpose:

This SOP describes the technique for qualitative and quantitative estimation of isoniazid in human plasma by Liquid Chromatography Tandem Mass Spectrometry (LCMS/MS).

2. Scope:

This SOP is limited to find out the concentration in  $\mu\text{g/ml}$  of isoniazid in human plasma by Liquid Chromatography Tandem Mass Spectrometry (LCMS/MS).

3. Responsibilities:

The head of the department is responsible for the medical care and welfare of all patients under her/his care. The task of performing estimation of isoniazid will be delegated to trained personnel who will perform this function.

4. Applicable rules, regulations and guidelines

- ICMR Good Clinical Laboratory Practices Guidelines 2021 (<http://icmr.nic.in/guidelines/GCLP.pdf>)

5. Reference to other applicable SOPs

- SOP No.24/01: Biomedical waste management.

Category : Patient Care  
Title : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).  
SOP No. : TDM12/01  
Date first effective: 1<sup>st</sup> January 2025 Review date: 31<sup>st</sup> December 2025  
Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

## 6. Detailed instructions:

### A. **Chemicals and Materials:**

- i. Pure Powder of Isoniazid
- ii. Pure powder of Phenacetin (IS)
- iii. Miscellaneous: : Methanol and Acetonitrile (HPLC Grade-Merck make), Formic Acid (HPLC Grade-Merck make), Ammonium acetate buffer (GR grade-Merck), Purified water.
- iv. LC column: C18, 150 mm x 4.6 mm, 5-micron particle size
- v. Glassware: Eppendorf vials (2.0mL), Glass Test Tubes (10ml), Conical Flasks (100, 250, 500, 1000 mL), Measuring Cylinders (100ml, 500ml), falcon tubes, volumetric flasks (10ml).

**B. Equipments:** High Performance Liquid Chromatography (Shimadzu), Mass Spectrometry (API 2000), Cooling Centrifuge 15,000 X g (Biofuge), Ultra sonicator machine (Imeco), Vortex machine (Spinage), Auto pipettes-eppendorf (10-100 $\mu$ L, 100-1000 $\mu$ L), Centrifuge (Remi-R23).

### C. **Preparation of solutions, mobile phase and calibration standards:**

#### a. **Solutions:**

- i. **Preparation of stock 5mM ammonium acetate buffer, pH 3.00:** - 0.3854 gm of ammonium acetate dissolved and diluted up to 1000 mL of purified water and adjusted pH 3.0 with formic acid solution.
- ii. **Preparation of aqueous methanol solution:** - 90 mL of methanol diluted up to 100 mL with 5mM ammonium acetate buffer solution and 0.1 mL of formic acid solution.
- iii. **. Preparation of 1: 1 Methanol: Purified water solution:**  
- 50 ml of Methanol diluted up to 100 mL with purified water.
- iv. **Preparation of 90% methanol solution (0.1 % formic acid solution):** 90 mL of methanol diluted up to 100 mL with purified water and add 100  $\mu$ L formic acid solution.

**Category** : Patient Care  
**Title** : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).  
SOP No. : TDM12/01

**Date first effective:** 1<sup>st</sup> January 2025      **Review date:** 31<sup>st</sup> December 2025  
**Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

- v. **Preparation of system suitability solution:** 975  $\mu$ L of mobile phase + 25  $\mu$ L of 1mg/mL of stock (1mg/mL)
- vi. **Mobile phase:** 5mM ammonium acetate buffer: Methanol: Acetonitrile (30:50:20)

**b. Calibration standards (Isoniazid):**

- i. **Preparation of aqueous stock standard solution (1 mg/ml) of isoniazid:** 10 mg of isoniazid hydrochloride dissolve and diluted up to 10 mL with 1:1 methanol: purified water solution and vortexes for 2 mins.
- ii. **Preparation of aqueous stock standard solution (100  $\mu$ g/ml) of isoniazid:** 1 mL of stock solution (1 mg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
- iii. **Preparation of aqueous standard solution (75  $\mu$ g/mL) of isoniazid:** 750  $\mu$ L of stock solution (1 mg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
- iv. **Preparation of aqueous standard solution (50  $\mu$ g/mL) of isoniazid:** 5 mL of stock solution (100  $\mu$ g/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
- v. **Preparation of aqueous standard solution (25  $\mu$ g/ml) of isoniazid:** 5 mL of stock solution (50  $\mu$ g/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)

**Category** : Patient Care  
**Title** : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).

SOP No. : TDM12/01

**Date first effective:** 1<sup>st</sup> January 2025 **Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

- vi. **Preparation of aqueous standard solution (10 µg/ml) of isoniazid:** 1 mL of stock solution (100 µg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
  - vii. **Preparation of aqueous standard solution of isoniazid (Tuning solution):** 1 mL of stock solution (10 µg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
  - viii. **Preparation of aqueous standard solution (0.5 µg/mL) of isoniazid:** 5 mL of stock solution (1 µg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
  - ix. **Preparation of aqueous standard solution (0.1 µg/mL) of isoniazid:** 1 mL of stock solution (10 µg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
  - x. **Preparation of aqueous standard solution (0.05 µg/mL) of isoniazid:** 5 mL of standard solution (0.1 µg/mL) diluted up to 10 mL with 90% methanol solution (0.1% formic acid)
- c. Preparation of Internal Standard (Phenacetin):**
- i. **Preparation of stock aqueous internal standard (1 mg/mL):** 10 mg of phenacetin dissolve and diluted up to 10 mL of methanol (100%) and vortexes for 2 mins.
  - ii. **Preparation of stock aqueous internal standard (100 µg/mL):** 1 mL of stock solution (1 mg/mL) diluted up to 10 mL with methanol (100%)

**Category** : Patient Care  
**Title** : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).  
SOP No. : TDM12/01

**Date first effective:** 1<sup>st</sup> January 2025 **Review date:** 31<sup>st</sup> December 2025  
**Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

- iii. **Preparation of working internal standard (10 µg/mL):** 1 mL of stock solution (100 µg/mL) diluted up to 10 mL with methanol (100%)
- iv. **Preparation of tuning solution (1 µg/mL):** 1 mL of stock solution (10 µg/mL) diluted up to 10 mL with 90% methanol (0.1% formic acid).

**d. Preparation of plasma Standard (Isoniazid):**

- i. **Preparation of stock plasma standard (100 µg/mL) of isoniazid:** 200 µL ml of aqueous stock solution (1 mg/mL) of isoniazid diluted with 1800 µL of blank plasma
- ii. **Preparation of plasma standard (75 µg/mL):** 750 µL of stock plasma standard (100 µg/mL) of isoniazid diluted with 250 µL of blank plasma.
- iii. **Preparation of plasma standard (50 µg/mL):** 500 µL of stock plasma standard (100 µg/mL) of isoniazid diluted with 500 µL mL of blank plasma.
- iv. **Preparation of plasma standard (25µg/mL):** 0.5 mL of stock plasma standard (50 µg/mL) of isoniazid diluted with 0.5 mL of blank plasma.
- v. **Preparation of plasma standard (10µg/ml):** 200 µL of stock plasma standard (100 µg/mL) of isoniazid diluted with 1800 µL of blank plasma.
- vi. **Preparation of plasma standard (5µg/ml):** 500 µL of plasma standard (10µg/mL) of isoniazid diluted with 500 µL of blank plasma.

Category : Patient Care  
Title : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).

SOP No. : TDM12/01

Date first effective: 1<sup>st</sup> January 2025

Review date: 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

- vii. **Preparation of plasma standard (1 µg/mL):** 200 µL of stock plasma standard (10µg/mL) of isoniazid diluted with 1800 µL of blank plasma.
  - viii. **Preparation of plasma standard (0.5µg/mL):** 500 µL of plasma standard (1 µg/mL) of isoniazid diluted with 500 µL of blank plasma.
  - ix. **Preparation of plasma standard (0.1µg/mL):** 200 µL of plasma standard (1.0 µg/mL) of isoniazid diluted with 1800 µL of blank plasma.
  - x. **Preparation of plasma standard (0.05µg/mL):** 500 µL of plasma standard (0.1 µg/mL) of isoniazid diluted with 500 µL of blank plasma
- e. **Preparation of quality control samples:**
- i. **Preparation of lower limit of quantification (LLOQ) sample (0.5µg/mL):** 500 µL of plasma standard (1 µg/mL) of isoniazid diluted with 500 µL of blank plasma.
  - ii. **Preparation of Middle Quality Control (MQC) sample (10 µg/ml):** 100 µL of stock plasma standard (100 µg/mL) of isoniazid diluted with 900 µL of blank plasma.
  - iii. **Preparation of High Quality Control (HQC) sample (40 µg/mL):** 400 µL of stock plasma standard (100 µg/mL) of isoniazid diluted with 600 µL of blank plasma.




**Category** : Patient Care  
**Title** : Procedure for estimation of isoniazid in Human Plasma  
by Liquid Chromatography Mass Spectrometry (LCMS).  
SOP No. : TDM12/01  
**Date first effective:** 1<sup>st</sup> January 2025 **Review date:** 31<sup>st</sup> December 2025  
Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

**f. Extraction procedure**

**Simple protein precipitation method**

100  $\mu$ L of BP/ plasma standard/QC + 50  $\mu$ L of (10  $\mu$ g/mL)  
working internal standard. Vortex for 1 minute.



**g. LCMS condition:**

**LC parameters**

- i. Injecting volume: 5 $\mu$ L
- ii. Flow rate: 0.350 mL/min
- iii. Column: C18, 150 mm x 4.6 mm, 5-micron particle size
- iv. Auto sampler temperature: 25° C
- v. Column oven temperature: 30° C
- vi. Total run time: 8 mins

**Category** : Patient Care  
**Title** : Procedure for estimation of isoniazid in Human Plasma  
 by Liquid Chromatography Mass Spectrometry (LCMS).  
 SOP No. : TDM12/01  
**Date first effective:** 1<sup>st</sup> January 2025 **Review date:** 31<sup>st</sup> December 2025  
**Department of Clinical Pharmacology, 1st Floor, New MS Building,**  
**Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

Comp name	Q1	Q3	Time (msec)	ID	DP	FP	EP	CE	CXP
INH	138.1	121.1	200	INH	20	310	5	21	5
Phenacetin	180.1	110.1	200	IS	23	250	8	30	5

**MS parameters:**

Scan type: MRM  
 Polarity: Positive  
 Curtain gas: 30 PSI  
 CAD: 10 PSI  
 Ion Spray voltage: 5500 v  
 Heater temperature: 425 °C  
 GS1: 50  
 GS2: 60

**1. Abbreviations:**

- i. **LCMS** = Liquid Chromatography Mass Spectrometry
- ii. **QC** = Quality Control
- iii. **LLOQ** : Lower Limit of Quantification
- iv. **MQC** : Middle Quality Control
- v. **HQC** : High Quality Control
- vi. **I.S.** = Internal Standard
- vii. **DP** = Declustering Potential
- viii. **FP** = Focusing Potential
- ix. **EP** = Entrance Potential
- x. **CE** = Collision Energy
- xi. **CXP** = Collision Exit Potential