

CURRENT

Category : Patient Care
Title : Procedure for Estimation of Hydroxyurea in Human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM-16/01

Date first effective: 1st January 2025 **Review date:** 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

Category : Patient Care

Title : Procedure for Estimation of Hydroxyurea in human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No.: TDM-16

Total pages: 6

Date first effective: 1st January 2025

Next Review date: 31st December 2025

Version: 01

Author:

Swati More
30/Dec/24
Swati More
Lab. Analyst

Dr. Bhaskar Krishnamurthy
Lab. In-charge

Signature with date:

Bsk
30/DEC/2024
Dr. Bhaskar Krishnamurthy
Assistant Professor,
Department of Clinical Pharmacology,
Seth GSMC and KEMH, Mumbai -400 012.

Reviewer:

Bell
30/Dec/2024
Dr. Mahesh Belhekar
Associate Professor
Department of Clinical Pharmacology
New MS Building, First Floor,
Seth GS Medical College and KEM Hospital
Acharya Donde Marg, Parel,
Mumbai - 400 012, India

Signature with date:

Approved by:

Dr. Nithya Gogtay
Professor and Head

Signature with date:

Dr. Nithya Gogtay
Professor & Head
Department of Clinical Pharmacology
1st Floor, MS Building,
Seth GS Medical College & KEM Hospit.
Parel, Mumbai - 400 012.

Category : Patient Care
Title : Procedure for Estimation of Hydroxyurea in Human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM-16/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
 Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

1. Table of Contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	Applicable rules/guidelines	3
5	Reference to other applicable SOPs	3
6	Detailed instructions	4
7	Abbreviations	6

Category : Patient Care
Title : Procedure for Estimation of Hydroxyurea in Human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM-16/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

1. Purpose:

This SOP describes the technique for qualitative and quantitative estimation of Hydroxyurea in human Plasma by High Performance Liquid Chromatography (HPLC).

2. Scope:

This SOP is limited to estimation of Hydroxyurea in human plasma by High Performance Liquid Chromatography (HPLC). **in the Department of Clinical Pharmacology, Seth GSMC and KEMH, Mumbai.**

3. Responsibilities:

The Head of the department is responsible for the medical care and welfare of all patients pertaining to TDM of Hydroxyurea. The task of performing estimation of Hydroxyurea will be delegated to trained personnel (laboratory technicians) who will perform this function.

4. Applicable rules, regulations and guidelines

- ICMR Good Clinical Laboratory Practices Guidelines 2008 (<http://icmr.nic.in/guidelines/GCLP.pdf>)

5. Reference to other applicable SOPs

- SOP No.24/02: Biomedical waste management.
- SOP No. TDM01/02: Collection and separation blood plasma for TDM
- SOP No. TDM05/02: Operation of High-Performance Liquid Chromatography

Category : Patient Care
Title : Procedure for Estimation of Hydroxyurea in Human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM-16/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

6. Detailed instructions

1. Preparation of standard and calibrator

a. Preparation of Stock Standard (Hydroxyurea 1mg/mL)

10 mg of pure powder of Hydroxyurea + 10 mL of Distilled water

b. Preparation of stock Internal standard (Carbamazepine 1mg/mL)

10mg of pure powder Carbamazepine + 10mL of Methanol

a. Working standard IS (25 ug/ml): 25 μ L of stock IS + 975 μ L of Methanol

c. Preparation of Stock (Xanthidrol 1mg/mL) :(10mg)

10 mg of pure powder of Xanthidrol + 10 mL Methanol

a. Working stock (100 ug/ml): 100 μ L of stock IS + 900 μ L of Methanol

d. Preparation of Stock Standard (Hydrochloric acid 1mg/mL) :(10mg)

10 mg of Hydrochloric acid + 10 mL Distilled Water

a. Working stock (100 ug/ml): 100 μ L of stock IS + 900 μ L of Distilled water

e. Preparation of Mobile Phase buffer (20mM):

Weigh 1.54165 gm of Ammonium acetate ($C_2H_7NO_2$) and dissolve in 1000 mL of distilled water

2. Preparation of plasma standards:

Concentration μ g/mL	20 (I)	10 (II)	5 (III)
Blank plasma	980 μ L	500 μ L	500 μ L
HU Stock 1mg/mL	20 μ L mix well	-----	---
		500 μ L stock(I)	500 μ L stock (II)

Category : Patient Care
Title : Procedure for Estimation of Hydroxyurea in Human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM-16/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

3. Extraction procedure

200 μ L of Plasma (blank plasma/standards/quality control/patients' sample) in 2mL Eppendorf tube.



Add 50 μ L of Internal standard (Working standard 25 ug/ml)



Add 600 μ L of Methanol



Add 100 μ L of Xanthinol (Working standard 100 ug/ml)



Add 50 μ L of HCl



The mixture and vortex 2 sec



Leave at room temperature protected from light for 5 min



Centrifuge the suspension at 8000 rpm for 10 mins



Directly Inject Supernant into HPLC for analysis.

4. Preparation of Mobile phase

- Take Ammonium acetate ($C_2H_7NO_2$) buffer (20mM): Acetonitrile; in the ratio of 700: 300 for 1000mL.
- Filter the mobile phase through 0.22-micron filter and sonicate for 15 minutes.

Category : Patient Care
Title : Procedure for Estimation of Hydroxyurea in Human Plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version: TDM-16/01

Date first effective: 1st January 2025

Review date: 31st December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

5. HPLC Conditions

- a. Injecting volume: 25 μ L
- b. Flow rate: 1.0mL/min.
- c. Wavelength: 240 nm (UV detector)
- d. Run Time: 25.00 min (approximately)
- e. Retention times for HU- 8.1- min, I.S 9.1 – min approximately.
- f. HPLC Column: Pico. Tag for Free Amino acid Analysis column (3.9 \times 300mm,).

6. Abbreviations:

- a. HPLC = High Performance Liquid Chromatography
- b. I.S. = Internal Standard
- c. HU = Hydroxyurea
- d. Std = Standard
- e. CBZ = Carbamazepine
- f. Xanth= Xanthidrol
- g. HCl = Hydrochloric acid