CURRENT

Category

: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by High Performance Liquid

Chromatography (HPLC) Eliminate of Nitrogen Gas SOP No. and Version: TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Category: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and

Carbamazepine (CBZ) in human blood plasma by High Performance Liquid

Chromatography (HPLC). Eliminate Nitrogen Gas

SOP No.: TDM 017

Total pages: 07

30/12/24

Date first effective: 1st January 2025

Next Review date: 31st December 2025

Version: 02

Author:

Lab. Analyst

Dr. Bhaskar Krishnamurthy

Lab. In charge

Signature with date:

631~ K 30/DEC/2024

Dr. Bhaskar Krishnamurthy Assistant Professor,

Department of Clinical Pharmacology, Seth GSMC and KEMH, Mumbai -400 012.

Reviewer:

Dr. Manesh Belhekar Associate Professor

De Mahesh N. Belhekar

Signature with date:

Associate Professor

Department of Clinical Pharmacology

New MS Building, First Floor, Soth GS Medical College and KEM Hospital Acharya Donde Marg, Parel,

Mumbai - 400 012. India

Approved by:

Dr. Nithya Gogtay **Professor and Head**

Signature with date: Professor & Head

Dr. Nithya Gogtay

Department of Clinical Pharmacology

1st Floor, MS Building,

Seth GS Medical College & KEM Hospital

Parel, Mumbai - 400 012.

Confidential

Page 1 of 6

CURRENT

Category

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by High Performance Liquid

Chromatography (HPLC) Eliminate of Nitrogen Gas

SOP No. and Version: TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

1. Table of Contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	Applicable rules/guidelines	3
5	Reference to other applicable SOPs	3
6	Detailed instructions	4
7	Abbreviations	7

Confidential

: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by High Performance Liquid Chromatography (HPLC) Eliminate of Nitrogen Gas

SOP No. and Version: TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

1. Purpose:

This SOP describes the technique for qualitative and quantitative estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by reverse phase High Performance Liquid Chromatography (HPLC).

2. Scope:

This SOP is limited to the estimation of PBT, PHT and CBZ in $\mu g/ml$ in human plasma by reverse phase High Performance Liquid Chromatography (HPLC).

3. Responsibilities:

The Head of the department is responsible for the medical care and welfare of all patients pertaining to TDM of PBT, PHT and CBZ. The task of performing estimation of PBT, PHT and CBZ will be delegated to trained personnel who will perform this function.

4. Applicable rules, regulations and guidelines

ICMR Good Clinical Laboratory Practices Guidelines 2008 (http://icmr.nic.in/guidelines/GCLP.pdf)

5. Reference to other applicable SOPs

- SOP No.24/02: Biomedical waste management.
- SOP No.TDM 01/02: Collection and separation of blood sample for TDM
- SOP No. TDM 05/02: Operation of High-Performance Liquid Chromatography

Confidential

: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and

Carbamazepine (CBZ) in human plasma by High Performance Liquid Chromatography (HPLC) Eliminate of Nitrogen Gas

SOP No. and Version: TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

6. Detailed instructions

- 1. Preparation of standards and calibrators
- 1.1. Preparation of Stock Standards
 - i. Preparation of Stock Standard (Phenobarbital-PBT, Phenytoin-PHT and Carbamazepine-CBZ) 1mg/ml:

(Phenobarbital-PBT) 1mg/ml	Add 10 mg of PBT in 10 ml of Methanol	Shake it well.
(Phenytoin-PHT) 1mg/ml	Add 10 mg of PHT in 10 ml of Methanol	Shake it well.
(Carbamazepine-CBZ) 1mg/ml	Add 10 mg of CBZ in 10 ml of Methanol	Shake it well.

iv. Preparation of Stock Internal Standard (10- Methoxy Carbamazepine I.S) 1mg/ml:

Add 10 mg of S.I.S in 10 ml of Methanol,	Shake it well.

V. Preparation of Working Internal Standard (5-ethyl-5-para-tolybarbituric acid-S.I.S) 1mg/ml:

Working standard (IS) I (100µg/ml)	10-Methoxy- Carbamazepine 1mg/ml 1mL	9 mL of Acetonitrile.	10ml
Working standard (IS) II (10μg/ml)	IS (I) 1mL	9 mL of Acetonitrile.	10ml

: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by High Performance Liquid Chromatography (HPLC) Eliminate of Nitrogen Gas

SOP No. and Version: TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

1.2. Preparation of Plasma Standards: 2ml Stock

Blank plasma	1750 μl	
(Phenobarbital-PBT)1mg/ml Stock	100 µl	
(Phenytoin-PHT) 1mg/ml	100 μΙ	
(Carbamazepine-CBZ) 1mg/ml	50 μl	

1.3. Preparation of Plasma Standards:

Standard I	0.9 ml blank plasma	0.1 ml of Stock
Standard II	0.8 ml blank plasma	0.2ml of Stock

1.4. Preparation of buffer for Mobile Phase:

- 0.2 M of KH_2PO_4 buffer at pH = 6.0
- Weigh 27.2 g of Potassium dihydrogen phosphate (KH_2PO_4) in 1000 mL of distilled water.
- Adjust pH = 6 with 1 M NaOH.

: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by High Performance Liquid

Chromatography (HPLC) Eliminate of Nitrogen Gas **SOP No. and Version:** TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

2. Extraction procedure

Add 100µl of plasma (blank/std/quality control/sample) in 2ml Plastic Eppendorf tube

Л

Add 100µl of internal standard (10µg/ml) to it Vortex Mix thoroughly for 30 sec

Л

Add 100µl of Methanol and Vortex Mix thoroughly for 2 min

 $\mathbf{1}$

Mix it well, white precipitate will be formed

Î

Extract the drug in organic phase by shaking it well (approx. 100 times)

Л

Centrifuge at 8000 rpm for 10 min

Û

Remove 100 μl of supernatant and add transfer into HPLC vials.

Inject into HPLC for analysis.

3. Preparation of Mobile phase

- a. Take 8.2 mL of buffer + 216.4 mL of Distilled Water + 185 mL of methanol + 90.6 mL of acetonitrile = 500 mL approximately for 40 samples.
- b. Filter the mobile phase through **0.22-micron** filter and sonicate for 15 minutes.

: Patient Care

Title

: Procedure for estimation of Phenobarbitone (PBT), Phenytoin (PHT) and Carbamazepine (CBZ) in human plasma by High Performance Liquid

Chromatography (HPLC) Eliminate of Nitrogen Gas

SOP No. and Version: TDM 17/02

Date first effective : 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

4. HPLC Conditions

a. Injecting volume: 10μL

b. Flow rate: 1.3 mL/min.

c. Wavelength: 255 nm

d. HPLC Column: C18 column

e. Run Time: 6.5 min (approximately)

f. Retention times for PBT- 3.4-4.0 min, I.S - 4.20 -4.50 min PHT- 4.6 - 5.0 and

CBZ - 5.3 - 6.0 min approximately.

5. Abbreviations:

- a. CBZ = Carbamazepine
- b. **HPLC** = High Performance Liquid Chromatography
- c. I.S. = Internal Standard
- d. **PBT** = Phenobarbital
- e. **PHT** = Phenytoin
- f. Std = Standard