CURRENT

Category

: Patient Care

Title

: Procedure for operation of High Performance Liquid Chromatography

(HPLC) machine.

SOP No. and Version: TDM 03/02

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Category

: Patient Care

Title

: Operation of High-Performance Liquid Chromatography (HPLC)

machine

SOP No.: TDM 03

Total pages: 05

Date first effective: 1st January 2025

Next Review date: 31st December 2025

Version: 02

Author:

Swati More

Lab. Analyst

Dr. Bhaskar Krishnamurthy Lab. In charge

> B1~1(30/0E(/2024

Dr. Bhaskar Krishnamurthy
Assistant Professor,
Department of Clinical Pharmacology,
Seth GSMC and KEMH, Mumbai -400 012.

Signature with date:

hely pedagh.

Reviewer:

Associate Professor

Mahesh N. Belhekar

Signature with date:

Associate Professor

Department of Clinical Pharmacology

Hew MS Building, First Floor,

Hew MS Building, First Floor,

Seth GS Medical College and KEM Hospital
Acharya Donde Marg, Parel,
Mumbai - 400 012. India

Approved by:

Dr. Nithya Gogtay Professor and Head

Signature with date

12 30112124

Dr. Nithya Gogtay

Department of Clinical Pharmacology

1st Floor, MS Building,

Seth GS Medical College & KEM Hospital,

Parel, Mumbai - 400 012.

Confidential

Page 1 of 5

CURRENT

Category

: Patient Care

Title

: Procedure for operation of High Performance Liquid Chromatography

(HPLC) machine.

SOP No. and Version: TDM 03/02

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

1. Table of Contents

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	Applicable rules/guidelines	3
5	Reference to other applicable SOPs	3
6	Detailed instructions	4
7	Abbreviations	5

Category

: Patient Care

Title

ころととととして

: Procedure for operation of High Performance Liquid Chromatography

(HPLC) machine.

SOP No. and Version: TDM 03/02

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

6. Detailed instructions

- Filter the solvents ((Distilled water, Acetonitrile: Distilled Water (70:30) or Methanol: Distilled on the rack in TDM assay lab) and mobile phase used for HPLC through 0.22-micron membrane filter
- 2. Sonicate all filtered solvents for 15 minutes in the sonicator
- 3. Switch on HPLC
- 4. Fix the column of the respective assay to HPLC (Direction of column should be in position from pump to detector).
- 5. Purge the HPLC pump with distilled water at 3 ml/min flow rate for 3 to 5 minutes to remove the air bubbles
- 6. Wash HPLC column with filtered and sonicated distilled water for 30 minutes at 1 ml/min flow rate.
- 7. Wash the HPLC column with filtered and sonicated Acetonitrile: Distilled Water (70:30) for VAL, VORI, VAN and LEV assay and Methanol: Distilled Water (70:30) for PBT, PHT, CBZ and LTG assay for 30 minutes at 1 ml/min flow rate
- 8. Prepare mobile phase as per the given in SOP of respective drug assay.
- 9. Filter the mobile phase and sonicate for 15 minutes
- 10. Saturate the column with mobile phase for 10 minutes at flow rate 1ml/min
- 11. Change the flow rate as per assay and keep for 10 minutes for stabilization
- 12. After stabilization commence with the injections according to assay protocol
- 13. After completion stop the pump and switch off the detector
- 14. Repeat steps 5 and 6
- 15. Switch off the HPLC

Category

: Patient Care

Title

: Procedure for operation of High Performance Liquid Chromatography

(HPLC) machine.

SOP No. and Version: TDM 03/02

Date first effective: 1st January 2025

Review date: 31st December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

1. Purpose:

This SOP describes the technique for the qualitative and quantitative estimation of drugs in human plasma using High Performance Liquid Chromatography (HPLC) machine.

2. Scope:

This SOP is limited to the technique used for Therapeutic drug monitoring (TDM) of drugs in human plasma by High Performance Liquid Chromatography (HPLC) machine.

3. Responsibilities:

The Lab technician and Lab In charge are responsible for the operation and maintenance of HPLC machine for TDM.

4. Applicable rules, regulations and guidelines

• ICMR Good Clinical Laboratory Practices Guidelines 2021 (http://icmr.nic.in/guidelines/GCLP.pdf)

5. Reference to other applicable SOPs

- SOP No.24/02: Biomedical waste management.
- SOP No.02/02: Estimation of PBT, PHT and CBZ
- SOP No. 03/02: Estimation of VAL
- SOP No. 04/02: Estimation of LTG
- SOP No. 07/02: Estimation of Li
- SOP No. 08/02: Estimation of LEV
- SOP No. 09/01: Estimation of PQ & CPQ
- SOP No. 10/01: Estimation of VAN
- SOP No. 11/01: Estimation of VORI

Category

: Patient Care

Title

: Procedure for operation of High Performance Liquid Chromatography

(HPLC) machine.

SOP No. and Version: TDM 03/02

Date first effective: 1st January 2025

Review date: 31st December 2025

Page 5 of 5

Department of Clinical Pharmacology, 1st Floor, New MS Building, Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

7. Abbreviations

- 1. HPLC-High Performance Liquid Chromatography
- 2. CBZ Carbamazepine
- 3. LTG Lamotrigine
- 4. PBT- Phenobarbitone
- 5. PHT- Phenytoin
- 6. VAL Valproate
- 7. LEV Levetiracetam
- 8. PQ-Primaquine
- 9. CPQ- Carboxyprimaquine
- 10. VAN- Vancomycin
- 11. VORI- Voriconazole