

# CURRENT

Category : Patient Care  
Title : Procedure for estimation of Valproate (VAL) in human plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version : TDM 04/02

Date first effective: 1<sup>st</sup> January 2025 Review date: 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

Category: Patient Care

Title: Procedure for estimation of Valproate (VAL) in human blood plasma by High Performance Liquid Chromatography (HPLC).

SOP No.: TDM 04/02

Total pages: 06

Date first effective: 1<sup>st</sup> January 2025

Next Review date: 31<sup>st</sup> December 2025

Version: 02

Author:

*Swati More*  
30/Dec/24  
Swati More  
Lab. Analyst

Dr. Bhaskar Krishnamurthy  
Lab In-charge

Signature with date:

*BK*  
30/DEC/2024

Dr. Bhaskar Krishnamurthy  
Assistant Professor,  
Department of Clinical Pharmacology,  
Seth GSMC and KEMH, Mumbai -400 012.

Reviewer:

*Mahesh*  
30/Dec/2024  
Dr. Mahesh Belhekar  
Associate Professor

Dr. Mahesh N. Belhekar  
Associate Professor

Signature with date:

Department of Clinical Pharmacology  
New MS Building, First Floor,  
Seth GS Medical College and KEM Hospital  
Acharya Donde Marg, Parel,  
Mumbai - 400 012. India

Approved by:

Dr. Nithya Gogtay  
Professor and Head

Signature with date:

*Nithya*  
30/12/24  
Dr. Nithya Gogtay  
Professor & Head  
Department of Clinical Pharmacology  
1<sup>st</sup> Floor, MS Building,  
Seth GS Medical College & KEM Hospital,  
Parel, Mumbai - 400 012.

**Category** : Patient Care  
**Title** : Procedure for estimation of Valproate (VAL) in human plasma by High Performance Liquid Chromatography (HPLC).  
**SOP No. and Version** : TDM 04/02

**Date first effective:** 1<sup>st</sup> January 2025  
**Review date:** 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

**1. Table of Contents**

No.	Contents	Page No.
1	Purpose	3
2	Scope	3
3	Responsibility	3
4	Applicable rules/guidelines	3
5	Reference to other applicable SOPs	3
6	Detailed instructions	4
7	Abbreviations	6

Category : Patient Care  
Title : Procedure for estimation of Valproate (VAL) in human plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version : TDM 04/02

Date first effective: 1<sup>st</sup> January 2025

Review date: 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

### 1. Purpose:

This SOP describes the technique for qualitative and quantitative estimation of Valproate in human plasma by High Performance Liquid Chromatography (HPLC).

### 2. Scope:

This SOP is limited to the estimation of Valproate in  $\mu\text{g/ml}$  in human plasma by High Performance Liquid Chromatography (HPLC).

### 3. Responsibilities:

The Head of the department is responsible for the medical care and welfare of all patients pertaining to TDM of Valproate. The task of performing estimation of Valproate will be delegated to trained personnel who will perform this function.

### 4. Applicable rules, regulations and guidelines

- ICMR Good Clinical Laboratory Practices Guidelines 2021 (<http://icmr.nic.in/guidelines/GCLP.pdf>)

### 5. Reference to other applicable SOPs

- SOP No.24/02: Biomedical waste management.
- SOP No. TDM01/02: Collection and separation of blood plasma for TDM
- SOP No. TDM05/02: Operation of High-Performance Liquid Chromatography

**Category** : Patient Care  
**Title** : Procedure for estimation of Valproate (VAL) in human plasma by High Performance Liquid Chromatography (HPLC).

**SOP No. and Version** : TDM 04/02

**Date first effective:** 1<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

**6. Detailed instructions**

**1. Preparation of standard and calibrator**

**a. Preparation of Stock Standard (Valproic acid 1mg/mL):(0.010g)**

10 mg of pure powder of valproate + 10 mL of distilled water

**b. Preparation of extraction buffer:**

0.952 g of potassium dihydrogen phosphate + 1.87 g of disodium hydrogen phosphate, dissolve in 50 mL of distilled water

**c. Preparation of extracting reagent:**

0.5g of 4-bromophenacyl bromide + 0.0125g of Dicyclohexane-18-crown-6, dissolve in 50mL of acetonitrile

**d. Preparation of stock Internal standard (IS): (Nonanoic acid)**

16.67 $\mu$ L (Nonanoic acid) of IS + 10mL of acetonitrile-stock IS

**e. Preparation of Working Internal standard (IS):**

Working standard I	1mL of stock IS	9 mL of Acetonitrile
Working standard II	1 mL of working standard I	9 mL of Acetonitrile
Working standard III	5 mL of working standard II	5 mL of Acetonitrile

Category : Patient Care  
Title : Procedure for estimation of Valproate (VAL) in human plasma by High Performance Liquid Chromatography (HPLC).

SOP No. and Version : TDM 04/02

Date first effective: 1<sup>st</sup> January 2025

Review date: 31<sup>st</sup> December 2025

Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.

e. Preparation of Mobile Phase buffer (0.003M):

a. Weigh 0.1021g of KH<sub>2</sub>PO<sub>4</sub> and dissolve in 250 mL of distilled water

Preparation of Valproic acid plasma standards:

	25µg/mL ( I)	50µg/mL (II)	100µg/mL(III)
Blank Plasma	375 µL	100µL	450 µL
valproic acid 1 mg/mL			50µL mix well transfer to
	125µL (100µg/mL)	100µL(100µg/mL)	

d. Extraction procedure

Add 0.1 mL of plasma (blank/std/quality control/sample)



0.025 mL of extracting buffer pipette in a clean dry 'U' shaped stoppered test tube



Add 0.25mL of internal standard to it, Mix thoroughly for 2 min.



Centrifuge at 2000 to 2500 rpm for 15mins



Transfer 200µL of supernatant in clean 'V' shaped test tube.

Add 50 µL of extracting solution



Esterify the content at 70°C for 15 minutes in water bath



Cool the tubes to room temperature.



And inject HPLC 50 µL

**Category** : Patient Care

**Title** : Procedure for estimation of Valproate (VAL) in human plasma by High Performance Liquid Chromatography (HPLC).

**SOP No. and Version** : TDM 04/02

**Date first effective:** 1<sup>st</sup> January 2025

**Review date:** 31<sup>st</sup> December 2025

**Department of Clinical Pharmacology, 1st Floor, New MS Building,  
Seth GS Medical College & KEM Hospital, Parel, Mumbai 400012.**

## **2. Preparation of Mobile phase**

- a. Take buffer: Acetonitrile in the ratio of 1:4
- b. Filter the mobile phase through 0.22-micron filter and sonicate for 15 minutes.

## **3. HPLC Conditions**

- a. Injecting volume: 50 $\mu$ L
- b. Flow rate: 1.4 mL/min.
- c. Wavelength: 254 nm
- d. Run Time: 9.00 min (approximately)
- e. Retention times for VAL- 5.0-6.0 min, I.S – 7.5 – 8.5 min approximately.
- f. HPLC Column: C18 column

## **4. Abbreviations:**

- a. HPLC = High Performance Liquid Chromatography
- b. I.S. = Internal Standard
- c. VAL = Valproate/Valproic Acid
- d. Std = Standard